
Reusable Pulse Oximetry (SpO2) Sensors Product Selection Document

Select type of the sensor from the below options and refer to *FMT Reusable Pulse Oximetry* (*SpO2*) *Sensors Compatibility Chart* to complete the exact part number.







RPS

Reusable – Pediatric Soft Finger Sensor Cable Length: 0.9 m or 3.0 m Application: Pediatric 10 - 40 kg Application Site: Index or other finger RIS Reusable – Infant Soft Finger Sensor Cable Length: 0.9 m or 3.0 m Application: Infant 3-15 kg Application Site: Index or other finger



THE REAL

RNS

Reusable – Neonate Wrap Soft Sensor Cable Length: 0.9 m or 3.0 m Application: Infant 3-15 kg; Neonate 1-4 kg Application Site: Infant/ great toe; Neonate/ foot

RWS

Reusable – Wrap Y Sensor Cable Length: 0.9 m or 3.0 m Application: Adult > 40 kg; Pediatric 10 - 40 kg; Infant 3-15 kg; Neonate 1-4 kg Application Site: Adult and Pediatric/ Index or other finger; Infant/ great toe; Neonate/ foot





RYS

Reusable –Y Sensor Cable Length: 0.9 m or 3.0 m Application: Adult > 40 kg; Pediatric 10 - 40 kg; Infant 3-15 kg; Neonate 1-4 kg Application Site: Adult and Pediatric/ Index or other finger; Infant/ great toe; Neonate/ foot

FMT ® is registered trademark of Metko Ltd. Other products, logos and company names mentioned herein maybe trademarks of their respective owners and are not to be taken as an endorsement or affiliation with FMT[®]. Specifications are subject to change without notice.

METKO is an ISO 9001:2000 and ISO13485:2003 certified manufacturer.